

EBP11653 - IMMUNIZING PEPTIDE

SKU: EBP11653



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002408.3; NP_001139438.1
Immunogen	Peptide with sequence ENESLQNGRKSTE-C, from the N Terminus of the protein sequence according to NP_002408.3; NP_001139438.1.
Product Comments	The peptide was used in the production of EBP11653 - Goat Anti-Antigen KI-67 / MKI67 (aa98-110) Antibody. It is currently untested in blocking applications.
Peptide Sequence	ENESLQNGRKSTE-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)