

EBP11654 - IMMUNIZING PEPTIDE

SKU: EBP11654



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_005564.1

Immunogen Peptide with sequence C-KESDLNGAQIKLRE, from the internal region of the protein sequence according to NP_005564.1.

Product Comments The peptide was used in the production of EBP11654 - Goat Anti-Lamin B1 (isoform 1) Antibody. It is currently untested in blocking applications.

Peptide Sequence C-KESDLNGAQIKLRE

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)