



EBP11693 - IMMUNIZING PEPTIDE

SKU: EBP11693



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002534.1; NP_001166103.1; NP_001166104.1
Immunogen	Peptide with sequence DLKIQTVKDQPDEK-C, from the N Terminus of the protein sequence according to NP_002534.1; NP_001166103.1; NP_001166104.1.
Product Comments	The peptide was used in the production of EBP11693 - Goat Anti-ORL1 / LOX1 / LOXIN Antibody. It is currently untested in blocking applications.
Peptide Sequence	DLKIQTVKDQPDEK-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)