

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11701 - IMMUNIZING PEPTIDE

SKU: EBP11701

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_000751.1; NP_724781.1; NP_758519.1 ID

Immunogen Peptide with sequence TQGDSILDCVPKDGQ, from the internal region of the protein sequence according to NP_000751.1; NP_724781.1; NP_758519.1.

The peptide was used in the production of EB11701 - Goat Anti-G-CSF receptor Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

TQGDSILDCVPKDGQ Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade