

EBP11702 - IMMUNIZING PEPTIDE

SKU: EBP11702



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage | |
| Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_000751.1; NP_724781.1; NP_758519.1 |
| Immunogen | Peptide with sequence C-PPELGKSPLTHYT, from the internal region of the protein sequence according to NP_000751.1; NP_724781.1; NP_758519.1. |
| Product | The peptide was used in the production of EBP11702 - Goat Anti-G-CSF receptor / CSF3R (aa551-563) Antibody. |
| Comments | It is currently untested in blocking applications. |
| Peptide Sequence | C-PPELGKSPLTHYT |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)