

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11702 - IMMUNIZING PEPTIDE**

**SKU:** EBP11702

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 000751.1; NP 724781.1; NP 758519.1 ID

Peptide with sequence C-PPELGKSPLTHYT, from the internal region of the protein sequence according to Immunogen

NP 000751.1; NP 724781.1; NP 758519.1.

The peptide was used in the production of EB11702 - Goat Anti-G-CSF receptor / CSF3R (aa551-563) Antibody. **Product** 

**Comments** It is currently untested in blocking applications.

**Peptide** 

C-PPELGKSPLTHYT Sequence

**Shipping** 

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade