

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11740 - IMMUNIZING PEPTIDE

SKU: EBP11740

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_000612.1; NP_001159419.1 ID

 $\label{eq:local_$

The peptide was used in the production of EB11740 - Goat Anti-5HT2A Receptor (aa296-309) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-QRSIHREPGSYTGR Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade