

EBP11747 - IMMUNIZING PEPTIDE

SKU: EBP11747



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001020537.2; NP_003367.4; NP_001020538.2; NP_001020539.2; etc.
Immunogen	Peptide with sequence C-RIKPHQGQHIGE, from the internal region of the protein sequence according to NP_001020537.2; NP_003367.4; NP_001020538.2; NP_001020539.2; etc..
Product	The peptide was used in the production of EB11747 - Goat Anti-VEGFA (aa288-299) Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-RIKPHQGQHIGE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)