

# EBP11748 - IMMUNIZING PEPTIDE

**SKU:** EBP11748



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_004439.2; NP_001005862.1; NP_001276867.1; NP_001276866.1; NP_001369713.1; NP_001369733.1
<b>Immunogen</b>	Peptide with sequence PREYVNARHCLP, from the internal region of the protein sequence according to NP_004439.2; NP_001005862.1; NP_001276867.1; NP_001276866.1; NP_001369713.1; NP_001369733.1.
<b>Product Comments</b>	The peptide was used in the production of EB11748 - Goat Anti-ERBB2 / HER2 Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	PREYVNARHCLP
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)