



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11748 - IMMUNIZING PEPTIDE

SKU: EBP11748

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_004439.2;\ NP_001005862.1;\ NP_001276867.1;\ NP_001276866.1;\ NP_001369713.1;\ NP_001369733.1$ ID

Peptide with sequence PREYVNARHCLP, from the internal region of the protein sequence according to NP_004439.2; Immunogen NP_001005862.1; NP_001276867.1; NP_001276866.1; NP_001369713.1; NP_001369733.1.

The peptide was used in the production of EB11748 - Goat Anti-ERBB2 / HER2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

PREYVNARHCLP Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade