

EBP11749 - IMMUNIZING PEPTIDE

SKU: EBP11749



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_004439.2; NP_001005862.1
Immunogen	Peptide with sequence TLIDTNRSRACHP, from the internal region of the protein sequence according to NP_004439.2; NP_001005862.1.
Product	The peptide was used in the production of EBP11749 - Goat Anti-ERBB2 / HER2 (aa182-194) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	TLIDTNRSRACHP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)