

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11752 - IMMUNIZING PEPTIDE

SKU: EBP11752

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001954.2; NP 001171601.1; NP 001171602.1 ID

Peptide with sequence SRQERVCNIEKNVS, from the internal region (near the N terminus) of the protein **Immunogen**

sequence according to NP 001954.2; NP 001171601.1; NP 001171602.1.

The peptide was used in the production of EB11752 - Goat Anti-Pro-EGF Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide SRQERVCNIEKNVS Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade