

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11757 - IMMUNIZING PEPTIDE

SKU: EBP11757

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_000563.1 ID

Immunogen Peptide with sequence QSENSCTHFPGNLPN, from the internal region of the protein sequence according to NP_000563.1.

The peptide was used in the production of EB11757 - Goat Anti-IL-10 (aa25-39) Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

QSENSCTHFPGNLPN

Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade