

# EBP11762 - IMMUNIZING PEPTIDE

**SKU:** EBP11762



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage** Shipped at ambient temperature, store at -20°C.

**Instructions**

**Immunogen** Peptide with sequence C-THVTYRKDTRYK, from the internal region of the protein sequence according to .

**Product** The peptide was used in the production of EBP11762 - Goat Anti-TIM-1 Antibody. It is currently untested in blocking applications.

**Peptide**  
**Sequence** C-THVTYRKDTRYK

**Shipping**  
**Instructions** Ambient

**Product**  
**Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)