

EBP11768 - IMMUNIZING PEPTIDE

SKU: EBP11768



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001103219.1; NP_766592.2
Immunogen	Peptide with sequence C-RRADGKPIARK, from the internal region of the protein sequence according to NP_001103219.1; NP_766592.2.
Product	The peptide was used in the production of EBP11768 - Goat Anti-Contactin 4 / Big-2 (mouse aa160-172)
Comments	Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RRADGKPIARK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)