

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11775 - IMMUNIZING PEPTIDE

SKU: EBP11775

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001158239.1 ID

 $\label{eq:local_$

The peptide was used in the production of EB11775 - Goat Anti-Robo3.1 (mouse aa1338-1351) Antibody. It is Product

Comments currently untested in blocking applications.

Peptide

C-PYGRPSFLSHGQGT Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade