

EBP11775 - IMMUNIZING PEPTIDE

SKU: EBP11775



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001158239.1
Immunogen	Peptide with sequence C-PYGRPSFLSHGQGT, from the C-terminus of the protein sequence according to NP_001158239.1.
Product Comments	The peptide was used in the production of EBP11775 - Goat Anti-Robo3.1 (mouse aa1338-1351) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PYGRPSFLSHGQGT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)