

EBP11776 - IMMUNIZING PEPTIDE

SKU: EBP11776



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | AAD11628.1 |
| Immunogen | Peptide with sequence KESQGRGRGLEACR, from the C-terminus of the protein sequence according to AAD11628.1. |
| Product Comments | The peptide was used in the production of EBP11776 - Goat Anti-Robo3.2 (mouse aa1303-1316) Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | KESQGRGRGLEACR |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)