

EBP11778 - IMMUNIZING PEPTIDE

SKU: EBP11778



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005714.2
Immunogen	Peptide with sequence SGKHYKGPEVSC, from the N Terminus of the protein sequence according to NP_005714.2.
Product Comments	The peptide was used in the production of EBP11778 - Goat Anti-TSPAN5 / NET-4 Antibody. It is currently untested in blocking applications.
Peptide Sequence	SGKHYKGPEVSC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)