

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11779 - IMMUNIZING PEPTIDE**

**SKU:** EBP11779

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 112189.2; NP 001121781.1 ID

Peptide with sequence C-SDIEAVKAGHH, from the C Terminus of the protein sequence according to Immunogen NP\_112189.2; NP\_001121781.1.

The peptide was used in the production of EB11779 - Goat Anti-TSPAN14 Antibody. It is currently untested in **Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Comments blocking applications.

**Peptide** 

C-SDIEAVKAGHH Sequence

Shipping

Grade

Instructions Ambient

**Product**