



EBP11812 - IMMUNIZING PEPTIDE

SKU: EBP11812



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006466.2; NP_001129406.1; NP_001129407.1; NP_001129408.1
Immunogen	Peptide with sequence C-TTQRYSDASKLREE, from the internal region of the protein sequence according to NP_006466.2; NP_001129406.1; NP_001129407.1; NP_001129408.1.
Product Comments	The peptide was used in the production of EBP11812 - Goat Anti-periostin (aa112-125) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TTQRYSDASKLREE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)