

EBP11822 - IMMUNIZING PEPTIDE

SKU: EBP11822



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006008.1; NP_001139319.1; NP_001139324.1; NP_001139323.1; NP_001139322.1; NP_001139321.1
Immunogen	Peptide with sequence HPSSETCNSIRLS, from the internal region of the protein sequence according to NP_006008.1; NP_001139319.1; NP_001139324.1; NP_001139323.1; NP_001139322.1; NP_001139321.1.
Product Comments	The peptide was used in the production of EBP11822 - Goat Anti-PROM1 / CD133 (aa310-322) Antibody. It is currently untested in blocking applications.
Peptide Sequence	HPSSETCNSIRLS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png