

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11822 - IMMUNIZING PEPTIDE

SKU: EBP11822

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

ID

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 006008.1; NP 001139319.1; NP 001139324.1; NP 001139323.1; NP 001139322.1; NP 001139321.1

Peptide with sequence HPSSETCNSIRLS, from the internal region of the protein sequence according to Immunogen

NP 006008.1; NP 001139319.1; NP 001139324.1; NP 001139323.1; NP 001139322.1; NP 001139321.1.

The peptide was used in the production of EB11822 - Goat Anti-PROM1 / CD133 (aa310-322) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

HPSSETCNSIRLS Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png