

EBP11847 - IMMUNIZING PEPTIDE

SKU: EBP11847



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_079765.3
Immunogen	Peptide with sequence C-QPTNTNEKTQLIRE, from the C Terminus of the protein sequence according to NP_079765.3.
Product Comments	The peptide was used in the production of EB11847 - Goat Anti-CSP-beta (mouse aa177-190) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QPTNTNEKTQLIRE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)