

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11871 - IMMUNIZING PEPTIDE

SKU: EBP11871

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_109652.3 ID

Immunogen Peptide with sequence C-RSLVQEGTGGKTQ, from the internal region of the protein sequence according to NP_109652.3.

The peptide was used in the production of EB11871 - Goat Anti-Prestin (aa399-411) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-RSLVQEGTGGKTQ Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade