

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11875-B - IMMUNIZING PEPTIDE

SKU: EBP11875-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 055018.2 ID

Peptide with sequence CELEPPAPESPMALLADPAPAAD., from the internal region of the protein sequence **Immunogen**

according to NP 055018.2.

The peptide was used in the production of EB11875-B - Goat Anti-AP2A1 (aa706-727), Biotinlyated Antibody. It **Product**

Comments is currently untested in blocking applications.

Peptide CELEPPAPESPMALLADPAPAAD.

Sequence **Shipping**

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade