

EBP11885 - IMMUNIZING PEPTIDE

SKU: EBP11885



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005808.3; NP_001157661.1; NP_001157666.1
Immunogen	Peptide with sequence C-TRGAVQSGVDKTKS, from the internal region of the protein sequence according to NP_005808.3; NP_001157661.1; NP_001157666.1.
Product	The peptide was used in the production of EBP11885 - Goat Anti-perilipin 3 / TIP47 (aa154-167) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-TRGAVQSGVDKTKS
Shipping	
Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)