

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11894 - IMMUNIZING PEPTIDE

SKU: EBP11894

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001013801.2 ID

Peptide with sequence C-DHITEEELKRFKY, from the internal region (near N Terminus) of the protein sequence

according to NP_001013801.2.

The peptide was used in the production of EB11894 - Goat Anti-Aim2 (rodent) Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-DHITEEELKRFKY Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade