



EBP11896 - IMMUNIZING PEPTIDE

SKU: EBP11896



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_079030.3
Immunogen	Peptide with sequence C-TLHITGQVTDSL D, from the internal region (near C terminus) of the protein sequence according to NP_079030.3.
Product Comments	The peptide was used in the production of EBP11896 - Goat Anti-PTCD2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TLHITGQVTDSL D
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)