

EBP11913 - IMMUNIZING PEPTIDE

SKU: EBP11913



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_997647.1; NP_002017.1; NP_997643.1; NP_997641.1; NP_997639.1; NP_473375.2

Immunogen Peptide with sequence KCERHTSVQTTSS, from the internal region (near N terminus) of the protein sequence according to NP_997647.1; NP_002017.1; NP_997643.1; NP_997641.1; NP_997639.1; NP_473375.2.

Product Comments The peptide was used in the production of EBP11913 - Goat Anti-fibronectin Antibody. It is currently untested in blocking applications.

Peptide Sequence KCERHTSVQTTSS

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png