

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11928 - IMMUNIZING PEPTIDE

SKU: EBP11928

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_003659.2; NP_620777.1 ID

Immunogen Peptide with sequence C-PQVLEAQRRHQKEKS, from the internal region of the protein sequence according to NP_003659.2; NP_620777.1.

The peptide was used in the production of EB11928 - Goat Anti-MK5 / MAPKAPK5 (aa192-206) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-PQVLEAQRRHQKEKS Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade