

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP11939 - IMMUNIZING PEPTIDE**

**SKU:** EBP11939

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 003274.3; NP 001119604.1; NP 001119605.1 ID

Peptide with sequence C-EKAQELSDWIHQLE, from the internal region of the protein sequence according to Immunogen

NP 003274.3; NP 001119604.1; NP 001119605.1.

The peptide was used in the production of EB11939 - Goat Anti-TNNT1 Antibody. It is currently untested in **Product** 

Comments blocking applications. **Peptide** 

C-EKAQELSDWIHQLE Sequence

Shipping

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade