

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP11975 - IMMUNIZING PEPTIDE

SKU: EBP11975

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_001899.1 ID

Product The peptide was used in the production of EB11975 - Goat Anti-Cathepsin B Antibody. It is currently untested in

Comments blocking applications.

Peptide

SDRICIHTNAHVSVE Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade