

EBP11976 - IMMUNIZING PEPTIDE

SKU: EBP11976



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001899.1; NP_001304166.1
Immunogen	Peptide with sequence SRPPCTGEGDTPK, from the internal region of the protein sequence according to NP_001899.1; NP_001304166.1.
Product	The peptide was used in the production of EBP11976 - Goat Anti-Cathepsin B (aa194-206) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	SRPPCTGEGDTPK
Shipping	
Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png