

# EBP11981 - IMMUNIZING PEPTIDE

SKU: EBP11981



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_004090.4; NP\_937837.1

**Immunogen** Peptide with sequence EKRHTRDSEAQR-C, from the N Terminus of the protein sequence according to NP\_004090.4; NP\_937837.1.

**Product Comments** The peptide was used in the production of EBP11981 - Goat Anti-Stomatin Antibody. It is currently untested in blocking applications.

**Peptide Sequence** EKRHTRDSEAQR-C

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)