

EBP11995 - IMMUNIZING PEPTIDE

SKU: EBP11995



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001186040.1; NP_005908.1; NP_001186041.1
Immunogen	Peptide with sequence C-PDVNHAKVKLQGK, from the internal region of the protein sequence according to NP_001186040.1; NP_005908.1; NP_001186041.1.
Product	The peptide was used in the production of EBP11995 - Goat Anti-MDH1 / MOR2 (aa211-223) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-PDVNHAKVKLQGK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)