

EBP12023 - IMMUNIZING PEPTIDE

SKU: EBP12023



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_073626.1; NP_758525.1
Immunogen	Peptide with sequence RPNRHPESCRASED, from the internal region of the protein sequence according to NP_073626.1; NP_758525.1.
Product Comments	The peptide was used in the production of EB12023 - Goat Anti-IL-25 / IL-17E (aa66-79) Antibody. It is currently untested in blocking applications.
Peptide Sequence	RPNRHPESCRASED
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)