

EBP12049 - IMMUNIZING PEPTIDE

SKU: EBP12049



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001106697.1; NP_001106699.1; NP_055732.2
Immunogen	Peptide with sequence C-SLYGQLDTTKLRD, from the internal region of the protein sequence according to NP_001106697.1; NP_001106699.1; NP_055732.2.
Product Comments	The peptide was used in the production of EBP12049 - Goat Anti-Netrin-G1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SLYGQLDTTKLRD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)