

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12060 - IMMUNIZING PEPTIDE

SKU: EBP12060

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_057297.2 ID

Immunogen Peptide with sequence C-DGQKSGDHLKRGR, from the internal region of the protein sequence according to NP_057297.2.

The peptide was used in the production of EB12060 - Goat Anti-TRIP4 / ASC1 (aa94-106) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-DGQKSGDHLKRGR Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade