

EBP12130 - IMMUNIZING PEPTIDE

SKU: EBP12130



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_065079.1
Immunogen	Peptide with sequence C-EDELDEETGDITQN, from the internal region of the protein sequence according to NP_065079.1.
Product Comments	The peptide was used in the production of EBP12130 - Goat Anti-VGLUT2 (aa521-534) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EDELDEETGDITQN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)