

EBP12156 - IMMUNIZING PEPTIDE

SKU: EBP12156



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001164185.1; NP_001164186.1; NP_001164187.1; NP_001164188.1; NP_001164189.1; NP_055382.2; NP_001164190.1; NP_001164191.1; NP_001164192.1
Immunogen	Peptide with sequence C-KQTSTFSKQTPHH, from the internal region of the protein sequence according to NP_001164185.1; NP_001164186.1; NP_001164187.1; NP_001164188.1; NP_001164189.1; NP_055382.2; NP_001164190.1; NP_001164191.1; NP_001164192.1.
Product Comments	The peptide was used in the production of EBP12156 - Goat Anti-BCAR1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KQTSTFSKQTPHH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)