

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12156 - IMMUNIZING PEPTIDE

SKU: EBP12156

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_001164185.1; NP_001164186.1; NP_001164187.1; NP_001164188.1; NP_001164189.1; NP_055382.2; Accession

NP_001164190.1; NP_001164191.1; NP_001164192.1 ID

Peptide with sequence C-KQTSTFSKQTPHH, from the internal region of the protein sequence according to

Immunogen NP 001164185.1; NP 001164186.1; NP 001164187.1; NP 001164188.1; NP 001164189.1; NP 055382.2;

NP 001164190.1; NP 001164191.1; NP 001164192.1.

The peptide was used in the production of EB12156 - Goat Anti-BCAR1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-KQTSTFSKQTPHH Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade