

EBP12229 - IMMUNIZING PEPTIDE

SKU: EBP12229



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001248.1; NP_001207417.1; NP_001207418.1, NP_001207419.1
Immunogen	Peptide with sequence C-NIRQQTPDKPSPN, from the internal region of the protein sequence according to NP_001248.1; NP_001207417.1; NP_001207418.1, NP_001207419.1.
Product Comments	The peptide was used in the production of EB12229 - Goat Anti-CDH13 / H cadherin (aa285-297) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NIRQQTPDKPSPN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)