

EBP12231 - IMMUNIZING PEPTIDE

SKU: EBP12231



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001014433.1; NP_057005.1; NP_001014840.1
Immunogen	Peptide with sequence TCPNEKVAKEIAR, from the internal region of the protein sequence according to NP_001014433.1; NP_057005.1; NP_001014840.1.
Product Comments	The peptide was used in the production of EBP12231 - Goat Anti-CutA protein (aa93-105) Antibody. It is currently untested in blocking applications.
Peptide Sequence	TCPNEKVAKEIAR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)