

EBP12255 - IMMUNIZING PEPTIDE

SKU: EBP12255



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000601.3; NP_001001389.1; NP_001001390.1; NP_001001391.1; NP_001001392.1; NP_001189484.1; NP_001189485.1; NP_001189486.1
Immunogen	Peptide with sequence C-NRDGTRYVQKGEYRT, from the internal region of the protein sequence according to NP_000601.3; NP_001001389.1; NP_001001390.1; NP_001001391.1; NP_001001392.1; NP_001189484.1; NP_001189485.1; NP_001189486.1.
Product Comments	The peptide was used in the production of EBP12255 - Goat Anti-CD44 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NRDGTRYVQKGEYRT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)