

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12255 - IMMUNIZING PEPTIDE

SKU: EBP12255

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

 $NP_000601.3; NP_001001389.1; NP_001001390.1; NP_001001391.1; NP_001001392.1; NP_001189484.1; NP_001001392.1; NP_001001001392.1; NP_001001392.1; NP_001001392.1; NP_001001392.1; NP_001001392.1; NP_001001392.1; NP_001001392$ Accession

NP_001189485.1; NP_001189486.1 ID

Peptide with sequence C-NRDGTRYVQKGEYRT, from the internal region of the protein sequence according to

Immunogen NP 000601.3; NP 001001389.1; NP 001001390.1; NP 001001391.1; NP 001001392.1; NP 001189484.1;

NP 001189485.1; NP 001189486.1.

The peptide was used in the production of EB12255 - Goat Anti-CD44 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-NRDGTRYVQKGEYRT Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade