

EBP12278 - IMMUNIZING PEPTIDE

SKU: EBP12278



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001138579.1; NP_006150.1; NP_001138581.1; NP_001138582.1
Immunogen	Peptide with sequence DPCQADTIRNDITKT, from the internal region (near C terminus) of the protein sequence according to NP_001138579.1; NP_006150.1; NP_001138581.1; NP_001138582.1.
Product Comments	The peptide was used in the production of EBP12278 - Goat Anti-NELL2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	DPCQADTIRNDITKT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)