

EBP12302 - IMMUNIZING PEPTIDE

SKU: EBP12302



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_055945.2

Immunogen Peptide with sequence C-DGERQFNLNGNSVPT, from the internal region of the protein sequence according to NP_055945.2.

Product Comments The peptide was used in the production of EB12302 - Goat Anti-TBC1D9 (aa444-458) Antibody. It is currently untested in blocking applications.

Peptide Sequence C-DGERQFNLNGNSVPT

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)