

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP12318 - IMMUNIZING PEPTIDE**

**SKU:** EBP12318

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_056150.1 ID

Immunogen Peptide with sequence C-SPKRSGKIPPKLHNH, from the internal region of the protein sequence according to NP\_056150.1.

The peptide was used in the production of EB12318 - Goat Anti-MED13L (aa361-375) Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

Peptide

C-SPKRSGKIPPKLHNH Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade