

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12320-B - IMMUNIZING PEPTIDE

SKU: EBP12320-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 001209.1; NP 996808.1; NP 001280571.1 ID

Peptide with sequence CSHLQHVKYKGQE., from the internal region of the protein sequence according to Immunogen

NP 001209.1; NP 996808.1; NP 001280571.1.

The peptide was used in the production of EB12320-B - Goat Anti-carbonic anhydrase XII (aa188-199), **Product**

Comments Biotinylated Antibody. It is currently untested in blocking applications.

Peptide

CSHLQHVKYKGQE. Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade